

DICOMRTTool_Tutorial

January 11, 2021

1 DICOM RT Tool Tutorial with Open-Access Data

This notebook demonstrates the various functions and utilities available in the Dicom RT tool Python package (https://github.com/brianmanderson/Dicom_RT_and_Images_to_Mask) by Anderson et. al. It serves as supplementary information for the Technical Paper titled: “Simple Python Module for Conversions between DICOM Images and Radiation Therapy Structures, Masks, and Prediction Arrays” . This notebook works through an example of publicly available brain tumor data of T1-w/FLAIR MRI sequences and corresponding RT structure files with multiple segmented regions of interest. Full information of the publicly available brain tumor data used in this notebook can be found at: https://figshare.com/articles/dataset/Data_from_An_Investigation_of_Machine_Learning_Methods_in_Delta_radiomics_Feature_Analysis/9943334. This notebook was written for easy accessibility for beginners to Python programming, medical imaging, and computational analysis. It should take no more than 10-15 minutes to run in it’s entirety from scratch. The notebook generates about 10 GB worth of files, so ensure you have adequate space to run it.

The notebook covers the following topics (click to go to section): 1. Getting the data 2. Reading in DICOM and RT struct files and converting to numpy array format 3. Saving arrays to nifti format and reloading them 4. Saving and loading numpy array files 5. Calculating radiomic features 6. Predictions To RT-Structure Example

The notebook assumes you have the following nested directory structure after running cells that download necessary data:

```
[ ]: """
Top-level directory/
  DICOMRTTool_manuscript.ipynb
  Example_Data/ <- Generated when you run the cells below
  /   Image_Data/
  /       Structure/ <- These correspond to the Pre-RT scans
  /           T1/
  /               Patient number/
  /                   RT Struc file (.dcm)
  /               T2FLAIR/
  /                   Patient number/
  /                       RT Struc file (.dcm)
  /           T1/
  /               Post1/
  /                   Patient number/
```

```

/          DICOM image files (.dcm)
/          Post2/
/          Patient number/
/          DICOM image files (.dcm)
/          Pre/
/          Patient number/
/          DICOM image files (.dcm) <- The images we care about
/          T2FLAIR/
/          Post1/
/          Patient number/
/          DICOM image files (.dcm)
/          Post2/
/          Patient number/
/          DICOM image files (.dcm)
/          Pre/
/          Patient number/
/          DICOM image files (.dcm) <- The images we care about
Data.zip <- Generated when you run the cells below, downloaded Figshare file
Nifti_Data/ <- Generated when you run the cells below
/          Image.nii
/          Mask.nii
/          MRN_Path_To_Iteration.xlsx
/          Overall_Data_Examples_(iteration)0.nii.gz
/          Overall_mask_Examples_y(iteration)0.nii.gz
Numpy_Data/ <- Generated when you run the cells below
/          image.npy
/          mask.npy
RT_Structures/ <- Generated when you run the cells below
/          RS_Test_UID.dcm
"""

```

```

[1]: %%capture
# Load or install the program, %%capture supresses print statements
!pip install DicomRTTool --upgrade
from DicomRTTool.ReaderWriter import DicomReaderWriter

```

```

[2]: # importing neccessary libraries

# file mangagment
import os
import zipfile
from six.moves import urllib

# array manipulation and plotting
import numpy as np
import pandas as pd
import matplotlib.pyplot as plt

```

```
# medical image manipulation
import SimpleITK as sitk
```

1.1 Part 1: Getting the data.

The RT struc files and their corresponding DICOM images can be in the same directory or different directories. Here we show a case where structure files and images are located in different directories. This is a good dataset to work with since its somewhat messy but coherent enough to show power of DICOMRTTool. Many files (pre-RT, post-RT at 2 timepoints) but only pre-RT T1 and FLAIR images have associated RT structure files. Downloading and unzipping the necessary files will take about 10 minutes on most CPUs and takes up about 8 GB of storage. One may visualize these DICOM images using a free commercially available DICOM viewer, such as Radiant (<https://www.radiantviewer.com/>).

```
[3]: %%time
data_path = os.path.join('.', 'Example_Data')
if not os.path.isdir(data_path): # create Example_data directory if it doesn't
    exist
    os.mkdir(data_path)

url_img = "https://ndownloader.figshare.com/files/20140100" # brain scans
filename_img = os.path.join(data_path, 'Data.zip')
if not os.path.exists(filename_img): # if zip file doesnt exist download
    print ("Retrieving zipped images...")
    print('Estimated download time is 5 minutes...')
    urllib.request.urlretrieve(url_img, filename_img)
    print('Finished downloading!')
else:
    print ("Zipped images already downloaded.")

if os.path.exists(filename_img): # If we downloaded the data
    if not os.path.exists(os.path.join(data_path, 'Image_Data')): # and it
    hasn't been unzipped
        print ("Unzipping images...")
        print('Estimated unzip time is 2 minutes!')
        z = zipfile.ZipFile(filename_img)
        z.extractall(data_path)
        print ("Done unzipping images.")

print("All required files downloaded and unzipped!") # print when done
```

```
Zipped images already downloaded.
All required files downloaded and unzipped!
Wall time: 1.13 ms
```

```
[4]: def display_slices(image, mask, skip=1):
    """
```

```

    Displays a series of slices in z-direction that contains the segmented
    ↪regions of interest.
    Ensures all contours are displayed in consistent and different colors.
    Parameters:
        image (array-like): Numpy array of image.
        mask (array-like): Numpy array of mask.
        skip (int): Only print every nth slice, i.e. if 3 only print every
    ↪3rd slice, default 1.
    Returns:
        None (series of in-line plots).
    """

    slice_locations = np.unique(np.where(mask != 0)[0]) # get indexes for where
    ↪there is a contour present
    slice_start = slice_locations[0] # first slice of contour
    slice_end = slice_locations[len(slice_locations)-1] # last slice of contour

    counter = 1

    for img_arr, contour_arr in zip(image[slice_start:slice_end+1],
    ↪mask[slice_start:slice_end+1]): # plot the slices with contours overlayed
    ↪ontop
        if counter % skip == 0: # if current slice is divisible by desired skip
    ↪amount
            masked_contour_arr = np.ma.masked_where(contour_arr == 0,
    ↪contour_arr)
            plt.imshow(img_arr, cmap='gray', interpolation='none')
            plt.imshow(masked_contour_arr, cmap='cool', interpolation='none',
    ↪alpha=0.5, vmin = 1, vmax = np.amax(mask)) # vmax is set as total number of
    ↪contours so same colors can be displayed for each slice
            plt.show()
            counter += 1

```

1.2 Part 2: Reading in DICOM and RT struct files and converting to numpy array format.

The principal on which this set of tools operates on is based on the DicomReaderWriter object. It is instantiated with the contours of interest (and associations) and can then be used to create numpy arrays of images and masks of the format [slices, width, height].

The following code logic is used to demonstrate searching a path and returning indices for matched structures and images (by UID) for arbitrary directory structures (DICOM image files and RT Struct files not in the same folder). If all necessary structure files are in the same folder as the corresponding images (by UID), one can alternatively use an os.walk through directories of interest and call DicomReaderWriter each time a folder is discovered. For example, I normally use a folder structure MRN -> date of image (pre,mid,post-RT) -> type of scan (MRI, CT, etc.) -> files (DICOM images + RT Struct). However, this approach calls the DicomReaderWriter iteratively,

which can be computationally taxing.

```
[5]: DICOM_path = os.path.join('.', 'Example_Data', 'Image_Data') # folder where  
      ↪ downloaded data was stored  
      print(DICOM_path)
```

```
.\Example_Data\Image_Data
```

This will walk through all of the folders, and using SimpleITK, will separate them based on SeriesInstanceUIDs.

```
[6]: %%time  
      Dicom_reader = DicomReaderWriter(description='Examples', arg_max=True)  
      print('Estimated 30 seconds, depending on number of cores present in your  
      ↪ computer')  
      Dicom_reader.walk_through_folders(DICOM_path) # need to define in order to use  
      ↪ all_roi method
```

```
Estimated 30 seconds, depending on number of cores present in your computer
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\006
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\001
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\004Loading from
```

```
.\Example_Data\Image_Data\Structure\T1\005
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\003
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\002
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\007
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\008
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\010Loading from
```

```
.\Example_Data\Image_Data\Structure\T1\009
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\011
```

```
Loading from .\Example_Data\Image_Data\Structure\T1\012
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\001
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\002
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\003Loading from
```

```
.\Example_Data\Image_Data\Structure\T2Flair\005Loading from
```

```
.\Example_Data\Image_Data\Structure\T2Flair\004
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\006Loading from
```

```
.\Example_Data\Image_Data\Structure\T2Flair\007
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\008
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\009Loading from
```

```
.\Example_Data\Image_Data\Structure\T2Flair\010
```

```
Loading from .\Example_Data\Image_Data\Structure\T2Flair\011Loading from
```

```
.\Example_Data\Image_Data\Structure\T2Flair\012
```

Loading from .\Example_Data\Image_Data\T1\002
Loading from .\Example_Data\Image_Data\T1\001
Loading from .\Example_Data\Image_Data\T1\003Loading from
.\Example_Data\Image_Data\T1\004Loading from .\Example_Data\Image_Data\T1\005

Loading from .\Example_Data\Image_Data\T1\006
Loading from .\Example_Data\Image_Data\T1\007
Loading from .\Example_Data\Image_Data\T1\008
Loading from .\Example_Data\Image_Data\T1\009
Loading from .\Example_Data\Image_Data\T1\010Loading from
.\Example_Data\Image_Data\T1\Post1\001

Loading from .\Example_Data\Image_Data\T1\Post1\002Loading from
.\Example_Data\Image_Data\T1\Post1\003

Loading from .\Example_Data\Image_Data\T1\Post1\004
Loading from .\Example_Data\Image_Data\T1\Post1\005
Loading from .\Example_Data\Image_Data\T1\Post1\006
Loading from .\Example_Data\Image_Data\T1\Post1\007
Loading from .\Example_Data\Image_Data\T1\Post1\008
Loading from .\Example_Data\Image_Data\T1\Post1\009
Loading from .\Example_Data\Image_Data\T1\Post1\010
Loading from .\Example_Data\Image_Data\T1\Post1\011
Loading from .\Example_Data\Image_Data\T1\Post1\012
Loading from .\Example_Data\Image_Data\T1\Post2\001
Loading from .\Example_Data\Image_Data\T1\Post2\002
Loading from .\Example_Data\Image_Data\T1\Post2\003
Loading from .\Example_Data\Image_Data\T1\Post2\004
Loading from .\Example_Data\Image_Data\T1\Post2\005
Loading from .\Example_Data\Image_Data\T1\Post2\006
Loading from .\Example_Data\Image_Data\T1\Post2\007
Loading from .\Example_Data\Image_Data\T1\Post2\008
Loading from .\Example_Data\Image_Data\T1\Post2\009
Loading from .\Example_Data\Image_Data\T1\Post2\010
Loading from .\Example_Data\Image_Data\T1\Post2\011
Loading from .\Example_Data\Image_Data\T1\Post2\012
Loading from .\Example_Data\Image_Data\T1\Pre\001
Loading from .\Example_Data\Image_Data\T1\Pre\002
Loading from .\Example_Data\Image_Data\T1\Pre\003
Loading from .\Example_Data\Image_Data\T1\Pre\004
Loading from .\Example_Data\Image_Data\T1\Pre\005
Loading from .\Example_Data\Image_Data\T1\Pre\006
Loading from .\Example_Data\Image_Data\T1\Pre\007
Loading from .\Example_Data\Image_Data\T1\Pre\008
Loading from .\Example_Data\Image_Data\T1\Pre\009
Loading from .\Example_Data\Image_Data\T1\Pre\010

Loading from .\Example_Data\Image_Data\T1\Pre\011
 Loading from .\Example_Data\Image_Data\T1\Pre\012
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\001
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\002
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\003
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\004
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\005
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\006
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\007
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\008
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\009
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\010
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\011
 Loading from .\Example_Data\Image_Data\T2Flair\Post1\012
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\001
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\002
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\003
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\004
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\005
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\006
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\007
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\008
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\009
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\010
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\011
 Loading from .\Example_Data\Image_Data\T2Flair\Post2\012
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\001
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\002
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\003
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\004
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\005
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\006
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\007
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\008
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\009
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\010
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\011
 Loading from .\Example_Data\Image_Data\T2Flair\Pre\012
 Compiling dictionaries together...
 Index 0, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Pre\002
 Index 1, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Post2\004
 Index 2, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Pre\006
 Index 3, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Post1\005
 Index 4, description Ax T2Flair Propeller at

.\Example_Data\Image_Data\T2Flair\Post1\006
Index 5, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Pre\007
Index 6, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\006
Index 7, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Pre\009
Index 8, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\010
Index 9, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\010
Index 10, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Pre\004
Index 11, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Pre\003
Index 12, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\002
Index 13, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\005
Index 14, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\009
Index 15, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\008
Index 16, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\009
Index 17, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\012
Index 18, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Pre\010
Index 19, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\011
Index 20, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\011
Index 21, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\004
Index 22, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\007
Index 23, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\007
Index 24, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\003
Index 25, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\012
Index 26, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\008
Index 27, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post1\001
Index 28, description Ax T2Flair Propeller at

.\Example_Data\Image_Data\T2Flair\Pre\005
 Index 29, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Pre\008
 Index 30, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Pre\001
 Index 31, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\003
 Index 32, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\001
 Index 33, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\002
 Index 34, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\005
 Index 35, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\005
 Index 36, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post1\008
 Index 37, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\006
 Index 38, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\006
 Index 39, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Pre\008
 Index 40, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\005
 Index 41, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\001
 Index 42, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Pre\009
 Index 43, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post1\009
 Index 44, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\007
 Index 45, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\006
 Index 46, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\007
 Index 47, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Post2\001
 Index 48, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\002
 Index 49, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Pre\004
 Index 50, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\003
 Index 51, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post2\004
 Index 52, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post2\001
 Index 53, description ax T1 3D 1MM fSPGR +C at
 .\Example_Data\Image_Data\T1\Post1\004
 Index 54, description Ax T2Flair Propeller at
 .\Example_Data\Image_Data\T2Flair\Pre\011
 Index 55, description Ax T2Flair Propeller at

```

.\Example_Data\Image_Data\T2Flair\Pre\012
Index 56, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\011
Index 57, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\008
Index 58, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Pre\010
Index 59, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\010
Index 60, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Pre\011
Index 61, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post1\011
Index 62, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post1\010
Index 63, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\007
Index 64, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post1\012
Index 65, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Pre\012
Index 66, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\009
Index 67, description ax T1 3D 1MM +c at .\Example_Data\Image_Data\T1\Post2\012
Index 68, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\002
Index 69, description ax T1 3D 1MM fSPGR +C at
.\Example_Data\Image_Data\T1\Post2\002
Index 70, description Ax T2Flair Propeller at
.\Example_Data\Image_Data\T2Flair\Post2\003
Index 71, description ax T1 3D 1MM fSPGR +C at
.\Example_Data\Image_Data\T1\Post2\003
Index 72, description None at None
Index 73, description None at None
74 unique series IDs were found. Default is index 0, to change use
set_index(index)
Wall time: 3.5 s

```

```

[7]: all_rois = Dicom_reader.return_rois(print_rois=True) # Return a list of all
      ↪ rois present, and print them

```

The following ROIs were found

```

rttempglioma
exprttempglioma
brainstem
dose 500[cgy]
dose 1000[cgy]
dose 1200[cgy]
gtvplus2
expltparrecgliom
ltparrecglioma
expltfrontrecao
ltfrontrecao
body
expltfrparrecgbm
ltfrparrecgbm
explttempglioma
lttempglioma
exprtfrontrecgbm

```

```

rtfrontrecgbm
expinftrttemprecg
inftrttempgbm
dose 2400[cgy]
expltfrontgbm
ltfrontgbm
exprttemprecglio
rttemprecglioma
rtfrontrecglioma
exprtfrontrecgli
brainstem1
eye, left
eye, right
chiasm
lens, left
lens, right
optic nerve, rig
optic nerve, lef
dose 2500[cgy]
exprttemprecgbm
rttemprecgbm
exprtfrparresxn
right_front_par_
abv
abv_roi

```

As we can see, these ROIs correspond to a variety of structures. In particular, we can see many GBM and glioma structures. Note GBM denotes glioblastoma multiforme (a high grade glioma).

```

[8]: # Print the locations of all RTs with a certain ROI name, automatically lower_
      ↪cased
      Dicom_reader.where_is_ROI(ROIName='BrAiNsTeM1')

```

Contours of brainstem1 are located:

```

.\Example_Data\Image_Data\Structure\T1\001\RS.CA1756_T13D.dcm
.\Example_Data\Image_Data\Structure\T1\011\RS.GF6065_T13D.dcm
.\Example_Data\Image_Data\Structure\T2Flair\001\RS.CA1756_T2Flair.dcm
.\Example_Data\Image_Data\Structure\T2Flair\011\RS.GF6065_T2Flairdcm.dcm
.\Example_Data\Image_Data\T1\001\RS.CA1756_T13D.dcm

```

```

[9]: Dicom_reader.which_indexes_have_all_rois() # Check to see which indexes have_
      ↪all of the rois we want
      # Since we haven't defined anything yet, it prompts you to input a list of_
      ↪contour names

```

You need to first define what ROIs you want, please use
`.set_contour_names_and_associations(roi_list)`

From these ROIs, we will look for those that describe the following regions of interest: tumor

(glioblastoma multiforme only) and high-dose area of radiation therapy.

```
[10]: Contour_Names = ['tumor', 'high_dose']
# Associations work as {'variant_name': 'desired_name'}
associations = {'dose 1000[cgy]': 'high_dose', 'dose 1200[cgy]': 'high_dose', #
↳new high dose ROI through association
                'exprtfrontrecgbm': 'tumor', 'rtfrontrecgbm': 'tumor',
↳'expltfrontgbm': 'tumor', # associating gbms to tumor
                'ltfrontgbm': 'tumor', 'infrttempgbm': 'tumor', 'rttemprecgbm':
↳'tumor',
                'exprttemprecgbm': 'tumor', 'expltfrparrecgbm': 'tumor',
↳'ltfrparrecgbm': 'tumor'}
```

```
[11]: Dicom_reader.set_contour_names_and_associations(Contour_Names=Contour_Names,
↳associations=associations)
```

```
Lacking ['tumor'] in index 0, location
.\Example_Data\Image_Data\T2Flair\Pre\002. Found ['rttempglioma',
'exprttempglioma', 'brainstem', 'dose 500[cgy]', 'dose 1000[cgy]', 'dose
1200[cgy]', 'gtvplus2']
Lacking ['tumor', 'high_dose'] in index 1, location
.\Example_Data\Image_Data\T2Flair\Post2\004. Found []
Lacking ['tumor'] in index 2, location
.\Example_Data\Image_Data\T2Flair\Pre\006. Found ['brainstem', 'dose 500[cgy]',
'dose 1200[cgy]', 'dose 1000[cgy]', 'expltparrecgliom', 'ltparrecglioma',
'gtvplus2']
Lacking ['tumor', 'high_dose'] in index 3, location
.\Example_Data\Image_Data\T2Flair\Post1\005. Found []
Lacking ['tumor', 'high_dose'] in index 4, location
.\Example_Data\Image_Data\T2Flair\Post1\006. Found []
Lacking ['tumor'] in index 5, location
.\Example_Data\Image_Data\T2Flair\Pre\007. Found ['brainstem', 'dose 500[cgy]',
'dose 1000[cgy]', 'dose 1200[cgy]', 'expltfrontrecao', 'ltfrontrecao', 'body',
'gtvplus2']
Lacking ['tumor', 'high_dose'] in index 6, location
.\Example_Data\Image_Data\T2Flair\Post2\006. Found []
Lacking ['tumor', 'high_dose'] in index 8, location
.\Example_Data\Image_Data\T2Flair\Post2\010. Found []
Lacking ['tumor', 'high_dose'] in index 9, location
.\Example_Data\Image_Data\T2Flair\Post1\010. Found []
Lacking ['tumor'] in index 10, location
.\Example_Data\Image_Data\T2Flair\Pre\004. Found ['brainstem', 'dose 1000[cgy]',
'dose 1200[cgy]', 'dose 500[cgy]', 'explttempglioma', 'lttempglioma',
'gtvplus2']
Lacking ['tumor', 'high_dose'] in index 12, location
.\Example_Data\Image_Data\T2Flair\Post1\002. Found []
Lacking ['tumor', 'high_dose'] in index 13, location
.\Example_Data\Image_Data\T2Flair\Post2\005. Found []
```

Lacking ['tumor', 'high_dose'] in index 14, location
 .\Example_Data\Image_Data\T2Flair\Post1\009. Found []
 Lacking ['tumor', 'high_dose'] in index 15, location
 .\Example_Data\Image_Data\T2Flair\Post2\008. Found []
 Lacking ['tumor', 'high_dose'] in index 16, location
 .\Example_Data\Image_Data\T2Flair\Post2\009. Found []
 Lacking ['tumor', 'high_dose'] in index 17, location
 .\Example_Data\Image_Data\T2Flair\Post2\012. Found []
 Lacking ['tumor', 'high_dose'] in index 19, location
 .\Example_Data\Image_Data\T2Flair\Post2\011. Found []
 Lacking ['tumor', 'high_dose'] in index 20, location
 .\Example_Data\Image_Data\T2Flair\Post1\011. Found []
 Lacking ['tumor', 'high_dose'] in index 21, location
 .\Example_Data\Image_Data\T2Flair\Post1\004. Found []
 Lacking ['tumor', 'high_dose'] in index 22, location
 .\Example_Data\Image_Data\T2Flair\Post1\007. Found []
 Lacking ['tumor', 'high_dose'] in index 23, location
 .\Example_Data\Image_Data\T2Flair\Post2\007. Found []
 Lacking ['tumor', 'high_dose'] in index 24, location
 .\Example_Data\Image_Data\T2Flair\Post1\003. Found []
 Lacking ['tumor', 'high_dose'] in index 25, location
 .\Example_Data\Image_Data\T2Flair\Post1\012. Found []
 Lacking ['tumor', 'high_dose'] in index 26, location
 .\Example_Data\Image_Data\T2Flair\Post1\008. Found []
 Lacking ['tumor', 'high_dose'] in index 27, location
 .\Example_Data\Image_Data\T2Flair\Post1\001. Found []
 Lacking ['tumor'] in index 29, location
 .\Example_Data\Image_Data\T2Flair\Pre\008. Found ['brainstem', 'dose 1000[cgy]',
 'dose 1200[cgy]', 'dose 500[cgy]', 'exprttempreglio', 'rttempreglioma',
 'gtvplus2']
 Lacking ['tumor'] in index 30, location
 .\Example_Data\Image_Data\T2Flair\Pre\001. Found ['rtfrontreglioma',
 'exprtfrontregli', 'brainstem1', 'dose 1000[cgy]', 'dose 1200[cgy]', 'dose
 500[cgy]', 'gtvplus2']
 Lacking ['tumor', 'high_dose'] in index 32, location
 .\Example_Data\Image_Data\T1\Post1\001. Found []
 Lacking ['tumor'] in index 33, location .\Example_Data\Image_Data\T1\Pre\002.
 Found ['exprttempglioma', 'rttempglioma', 'brainstem', 'dose 500[cgy]', 'dose
 1000[cgy]', 'dose 1200[cgy]', 'body', 'gtvplus2']
 Lacking ['tumor', 'high_dose'] in index 34, location
 .\Example_Data\Image_Data\T1\Post2\005. Found []
 Lacking ['tumor', 'high_dose'] in index 36, location
 .\Example_Data\Image_Data\T1\Post1\008. Found []
 Lacking ['tumor', 'high_dose'] in index 37, location
 .\Example_Data\Image_Data\T1\Post1\006. Found []
 Lacking ['tumor', 'high_dose'] in index 38, location
 .\Example_Data\Image_Data\T1\Post2\006. Found []
 Lacking ['tumor'] in index 39, location .\Example_Data\Image_Data\T1\Pre\008.

Found ['eye, left', 'eye, right', 'brainstem', 'chiasm', 'lens, left', 'lens, right', 'rttempreglioma', 'expttempreglio', 'dose 500[cgy]', 'dose 1000[cgy]', 'dose 1200[cgy]', 'gtvplus2', 'body']
 Lacking ['tumor', 'high_dose'] in index 40, location
 .\Example_Data\Image_Data\T1\Post1\005. Found []
 Lacking ['tumor'] in index 41, location .\Example_Data\Image_Data\T1\Pre\001.
 Found ['exptfrontregli', 'rtfrontreglioma', 'brainstem1', 'dose 500[cgy]', 'dose 1200[cgy]', 'dose 1000[cgy]', 'body', 'gtvplus2']
 Lacking ['tumor', 'high_dose'] in index 42, location
 .\Example_Data\Image_Data\T1\Pre\009. Found []
 Lacking ['tumor', 'high_dose'] in index 43, location
 .\Example_Data\Image_Data\T1\Post1\009. Found []
 Lacking ['tumor', 'high_dose'] in index 44, location
 .\Example_Data\Image_Data\T1\Post1\007. Found []
 Lacking ['tumor'] in index 45, location .\Example_Data\Image_Data\T1\Pre\006.
 Found ['eye, left', 'lens, right', 'lens, left', 'ltparreglioma', 'chiasm', 'optic nerve, rig', 'optic nerve, lef', 'expltparregliom', 'brainstem', 'eye, right', 'dose 500[cgy]', 'dose 1000[cgy]', 'dose 1200[cgy]', 'gtvplus2']
 Lacking ['tumor'] in index 46, location .\Example_Data\Image_Data\T1\Pre\007.
 Found ['brainstem', 'chiasm', 'expltfrontrecao', 'eye, left', 'eye, right', 'lens, left', 'lens, right', 'ltfrontrecao', 'optic nerve, lef', 'optic nerve, rig', 'dose 2500[cgy]', 'dose 500[cgy]', 'dose 1000[cgy]', 'dose 1200[cgy]', 'gtvplus2', 'body']
 Lacking ['tumor', 'high_dose'] in index 47, location
 .\Example_Data\Image_Data\T2Flair\Post2\001. Found []
 Lacking ['tumor', 'high_dose'] in index 48, location
 .\Example_Data\Image_Data\T1\Post1\002. Found []
 Lacking ['tumor'] in index 49, location .\Example_Data\Image_Data\T1\Pre\004.
 Found ['brainstem', 'dose 1000[cgy]', 'dose 1200[cgy]', 'dose 500[cgy]', 'explttempglioma', 'lttempglioma', 'gtvplus2']
 Lacking ['tumor', 'high_dose'] in index 50, location
 .\Example_Data\Image_Data\T1\Post1\003. Found []
 Lacking ['tumor', 'high_dose'] in index 51, location
 .\Example_Data\Image_Data\T1\Post2\004. Found []
 Lacking ['tumor', 'high_dose'] in index 52, location
 .\Example_Data\Image_Data\T1\Post2\001. Found []
 Lacking ['tumor', 'high_dose'] in index 53, location
 .\Example_Data\Image_Data\T1\Post1\004. Found []
 Lacking ['tumor'] in index 55, location
 .\Example_Data\Image_Data\T2Flair\Pre\012. Found ['brainstem', 'chiasm', 'dose 1000[cgy]', 'dose 1200[cgy]', 'dose 500[cgy]', 'exptfrparresxn', 'gtvplus2', 'right_front_par_']
 Lacking ['tumor', 'high_dose'] in index 56, location
 .\Example_Data\Image_Data\T1\Post2\011. Found []
 Lacking ['tumor', 'high_dose'] in index 57, location
 .\Example_Data\Image_Data\T1\Post2\008. Found []
 Lacking ['tumor', 'high_dose'] in index 59, location
 .\Example_Data\Image_Data\T1\Post2\010. Found []

```

Lacking ['tumor', 'high_dose'] in index 61, location
.\Example_Data\Image_Data\T1\Post1\011. Found []
Lacking ['tumor', 'high_dose'] in index 62, location
.\Example_Data\Image_Data\T1\Post1\010. Found []
Lacking ['tumor', 'high_dose'] in index 63, location
.\Example_Data\Image_Data\T1\Post2\007. Found []
Lacking ['tumor', 'high_dose'] in index 64, location
.\Example_Data\Image_Data\T1\Post1\012. Found []
Lacking ['tumor'] in index 65, location .\Example_Data\Image_Data\T1\Pre\012.
Found ['brainstem', 'chiasm', 'exptfrparresxn', 'right_front_par_', 'body',
'dose 500[cgy]', 'dose 1000[cgy]', 'dose 1200[cgy]', 'abv', 'abv_roi',
'gtvplus2']
Lacking ['tumor', 'high_dose'] in index 66, location
.\Example_Data\Image_Data\T1\Post2\009. Found []
Lacking ['tumor', 'high_dose'] in index 67, location
.\Example_Data\Image_Data\T1\Post2\012. Found []
Lacking ['tumor', 'high_dose'] in index 68, location
.\Example_Data\Image_Data\T2Flair\Post2\002. Found []
Lacking ['tumor', 'high_dose'] in index 69, location
.\Example_Data\Image_Data\T1\Post2\002. Found []
Lacking ['tumor', 'high_dose'] in index 70, location
.\Example_Data\Image_Data\T2Flair\Post2\003. Found []
Lacking ['tumor', 'high_dose'] in index 71, location
.\Example_Data\Image_Data\T1\Post2\003. Found []

```

Note: The module is printing “Found []” because many of the scans (post-1 and post-2 RT) do not have associated structure files. The module recognizes these images exist (unique UUIDs) but associated structure files cannot be located for them.

```

[12]: indexes = Dicom_reader.which_indexes_have_all_rois() # Check to see which
      ↪ indexes have all of the rois we want, now we can see indexes

```

The following indexes have all ROIs present

```

Index 7, located at .\Example_Data\Image_Data\T2Flair\Pre\009
Index 11, located at .\Example_Data\Image_Data\T2Flair\Pre\003
Index 18, located at .\Example_Data\Image_Data\T2Flair\Pre\010
Index 28, located at .\Example_Data\Image_Data\T2Flair\Pre\005
Index 31, located at .\Example_Data\Image_Data\T1\Pre\003
Index 35, located at .\Example_Data\Image_Data\T1\Pre\005
Index 54, located at .\Example_Data\Image_Data\T2Flair\Pre\011
Index 58, located at .\Example_Data\Image_Data\T1\Pre\010
Index 60, located at .\Example_Data\Image_Data\T1\Pre\011

```

```


[13]: pt_idx = indexes[-1]
      Dicom_reader.set_index(pt_idx) # This index has all the structures,
      ↪ corresponds to pre-RT T1-w image for patient 011
      Dicom_reader.get_images_and_mask() # Load up the images and mask for the
      ↪ requested index



```



```
Loading images for ax T1 3D 1MM +c at  
.\Example_Data\Image_Data\T1\Pre\011
```



```
[14]: image = Dicom_reader.ArrayDicom # image array  
mask = Dicom_reader.mask # mask array  
dicom_sitk_handle = Dicom_reader.dicom_handle # SimpleITK image handle  
mask_sitk_handle = Dicom_reader.annotation_handle # SimpleITK mask handle
```



```
[15]: n_slices_skip = 4  
display_slices(image, mask, skip = n_slices_skip) # visualize that our  
↪ segmentations were succesfully convereted
```











Note: Cyan color denotes tumor while magenta denotes surrounding area of high-dose radiation. Only displaying 7 slices.

1.3 Part 3: Saving arrays to nifti format.

If you want to use a manual approach, you can view the nifti files easily after running `get_images_and_mask()`. Saving files as nifti is advisable since spacing information is preserved.

```
[16]: nifti_path = os.path.join('.', 'Example_Data', 'Nifti_Data') # nifti subfolder
      if not os.path.exists(nifti_path):
          os.makedirs(nifti_path)
```

```
[17]: dicom_sitk_handle = Dicom_reader.dicom_handle # SimpleITK image handle
      mask_sitk_handle = Dicom_reader.annotation_handle # SimpleITK mask handle
      sitk.WriteImage(dicom_sitk_handle, os.path.join(nifti_path, 'Image.nii'))
      sitk.WriteImage(mask_sitk_handle, os.path.join(nifti_path, 'Mask.nii'))
```

One can also use the built in `.write_parallel` attribute to generate nifti files for all relevant pairs the `DicomReaderWriter` object has found/generated. In this case there are 9 image/mask pairs for unique UIDs that contain all contours we are interested in. Note a corresponding log excel file in the specified output path. The nifti files are written in the following format: “Overall_Data_{description}_ {iteration}.nii.gz” (image) or “Overall_mask_{description}_ y{iteration}.nii.gz” (mask).



```
[18]: %%time
      %%capture
      Dicom_reader.write_parallel(out_path = nifti_path, excel_file = os.path.
      ↪join(nifti_path, '.', 'MRN_Path_To_Iteration.xlsx'))
```



Wall time: 305 ms



We can now reload the nifti files and display them to check that nothing went wrong. You can inspect the other converted files by changing the numerical suffix as per the excel log file ('MRN_Path_To_Iteration.xlsx').



```
[19]: nifti_image = sitk.ReadImage(os.path.join(nifti_path, "Overall_Data_Examples_8.
      ↪nii.gz")) # reload image
      image = sitk.GetArrayFromImage(nifti_image)
      nifti_mask = sitk.ReadImage(os.path.join(nifti_path, "Overall_mask_Examples_y8.
      ↪nii.gz")) # reload mask
      mask = sitk.GetArrayFromImage(nifti_mask)
```


```
[20]: display_slices(image, mask, skip = n_slices_skip) # visualize that our
      ↪segmentations were succesfully converted from nifti
```











1.4 Part 4: Saving and loading numpy files for later use.

Finally we can save the numpy arrays themselves to files for later use (so you don't have to re-instantiate the computationally expensive DicomReaderWriter object) and subsequently re-load the numpy arrays.

```
[21]: numpy_path = os.path.join(data_path, 'Numpy_Data') # go into numpy subfolder
      if not os.path.exists(numpy_path):
          os.makedirs(numpy_path)
```

```
[22]: np.save(os.path.join(numpy_path, 'image'), image) # save the arrays
      np.save(os.path.join(numpy_path, 'mask'), mask)
```

```
[23]: image = np.load(os.path.join(numpy_path, 'image.npy')) # load the arrays
      mask = np.load(os.path.join(numpy_path, 'mask.npy'))
```

1.5 Part 5: Radiomics Use-case Example.

Here we use the popular open-source radiomics library PyRadiomics (<https://pyradiomics.readthedocs.io/en/latest/>) to calculate radiomic features for our ROIs. In this case, we only calculate a limited number features from the tumor as an illustrative example.

```
[24]: try:
      from radiomics import featureextractor
    except:
        !pip install pyradiomics
```

```
from radiomics import featureextractor
```

```
[25]: pd.set_option('display.max_columns', None) # show all columns
```

```
[26]: %%time
# note: need sitk images (sitk.ReadImage(nifti file)) to plug into PyRadiomics,
↳ preserves spacing

ROI_index = 1 # index for tumor
nifti_mask_tumor = sitk.BinaryThreshold(nifti_mask, lowerThreshold=ROI_index,
↳ upperThreshold=ROI_index) # select only ROI of interest

params = {} # can edit in more params as necessary
extractor = featureextractor.RadiomicsFeatureExtractor(**params) # instantiate
↳ extractor with parameters
extractor.disableAllFeatures() # in case where only want some features, can
↳ delete disable/enable lines if you want default
extractor.enableFeatureClassByName('firstorder')
extractor.enableFeatureClassByName('glcm')
features = {} # empty dictionary
features = extractor.execute(nifti_image, nifti_mask_tumor) # unpack results
↳ into features dictionary
df = pd.DataFrame({k: [v] for k, v in features.items()}) # put dictionary into
↳ a dataframe
```

GLCM is symmetrical, therefore Sum Average = 2 * Joint Average, only 1 needs to be calculated

Wall time: 4.06 s

```
[27]: df # display dataframe to inspect features
```

```
[27]: diagnostics_Versions_PyRadiomics diagnostics_Versions_Numpy \
0 v3.0.1 1.19.5

diagnostics_Versions_SimpleITK diagnostics_Versions_PyWavelet \
0 2.0.2 1.1.1

diagnostics_Versions_Python \
0 3.6.8

diagnostics_Configuration_Settings \
0 {'minimumROIDimensions': 2, 'minimumROISize': ...

diagnostics_Configuration_EnabledImageTypes \
0 {'Original': {}}

diagnostics_Image-original_Hash \
```

```

0  adcd47a617bf0a6906a361ba78d7e86388bc2fc9

    diagnostics_Image-original_Dimensionality \
0                                     3D

    diagnostics_Image-original_Spacing \
0  (0.5859000086784363, 0.5859000086784363, 1.0)

    diagnostics_Image-original_Size  diagnostics_Image-original_Mean \
0          (512, 512, 192)                                138.125106

    diagnostics_Image-original_Minimum  diagnostics_Image-original_Maximum \
0                                     0.0                                1936.0

    diagnostics_Mask-original_Hash \
0  ca99975f7f6ce5b9da272d5b065bb0fd8d40810e

    diagnostics_Mask-original_Spacing \
0  (0.5859000086784363, 0.5859000086784363, 1.0)

    diagnostics_Mask-original_Size  diagnostics_Mask-original_BoundingBox \
0          (512, 512, 192)                                (167, 225, 68, 37, 50, 16)

    diagnostics_Mask-original_VoxelNum  diagnostics_Mask-original_VolumeNum \
0                                     10905                                1

    diagnostics_Mask-original_CenterOfMassIndex \
0  (184.79257221458047, 254.31728564878497, 75.06...

    diagnostics_Mask-original_CenterOfMass \
0  (-41.43703106818862, 1.4954971831525938, -7.49...

    original_firstorder_10Percentile  original_firstorder_90Percentile \
0                                     659.0                                1042.0

    original_firstorder_Energy  original_firstorder_Entropy \
0          8464051347.0                                4.503961941913401

    original_firstorder_InterquartileRange  original_firstorder_Kurtosis \
0                                     188.0                                3.056529328502632

    original_firstorder_Maximum  original_firstorder_MeanAbsoluteDeviation \
0          1292.0                                114.14490366028481

    original_firstorder_Mean  original_firstorder_Median \
0          869.2987620357634                                888.0

```

```

original_firstorder_Minimum original_firstorder_Range \
0 311.0 981.0

original_firstorder_RobustMeanAbsoluteDeviation \
0 80.00459220139915

original_firstorder_RootMeanSquared original_firstorder_Skewness \
0 881.000814067669 -0.5230388203250016

original_firstorder_TotalEnergy original_firstorder_Uniformity \
0 2905529560.251311 0.051218440447186646

original_firstorder_Variance original_glcm_Autocorrelation \
0 20482.096710984642 579.312440764934

original_glcm_ClusterProminence original_glcm_ClusterShade \
0 42908.67954319103 -693.8616497595558

original_glcm_ClusterTendency original_glcm_Contrast \
0 118.47853188255739 8.74050977809287

original_glcm_Correlation original_glcm_DifferenceAverage \
0 0.8617988560957099 2.128732075933059

original_glcm_DifferenceEntropy original_glcm_DifferenceVariance \
0 2.6510141072718554 3.9429016300371384

original_glcm_Id original_glcm_Idm original_glcm_Idmn \
0 0.4588122607941779 0.3890481464071379 0.9946965137645681

original_glcm_Idn original_glcm_Imc1 original_glcm_Imc2 \
0 0.9515499470890418 -0.25514812891031574 0.9358524694499407

original_glcm_InverseVariance original_glcm_JointAverage \
0 0.37803667377261313 23.491920944502194

original_glcm_JointEnergy original_glcm_JointEntropy original_glcm_MCC \
0 0.006641162131651515 7.807733854014718 0.8696978190928191

original_glcm_MaximumProbability original_glcm_SumAverage \
0 0.01625971801846341 46.98384188900442

original_glcm_SumEntropy original_glcm_SumSquares
0 5.42170370906024 31.80476041516255

```

Numerical results for radiomic features shown here are consistent with importing nifti files as image and label map in 3D Slicer (<https://www.slicer.org/>) and using Radiomics extension (<https://www.slicer.org/wiki/Documentation/Nightly/Extensions/Radiomics>).

1.6 Part 6: Predictions To RT-Structure Example

Here we will provide a simple example for converting a predicted NumPy array of a square into a Dicom RT-Structure file

```
[28]: RT_path = os.path.join('Example_Data', 'RT_Structures')
      if not os.path.exists(RT_path):
          os.makedirs(RT_path)
```

First, we will create a fake prediction, it will be the same size as the image NumPy array

```
[29]: image = Dicom_reader.ArrayDicom
```

Now, deep learning model typically create segmentations in the format of (z_images, rows, cols, # of classes)

```
[30]: def create_circular_mask(h, w, center=None, radius=None):



      if center is None: # use the middle of the image
          center = (int(w/2), int(h/2))
      if radius is None: # use the smallest distance between the center and image
          ↪walls
          radius = min(center[0], center[1], w-center[0], h-center[1])



      Y, X = np.ogrid[:h, :w]
      dist_from_center = np.sqrt((X - center[0])**2 + (Y-center[1])**2)


      mask = dist_from_center <= radius
      return mask
```

```
[31]: predictions = np.zeros(image.shape + (4,)) # Four classes: background, square,
      ↪circle, target
      predictions.shape
      predictions[75:80, 250:350, 100:200, 1] = 1 # Here we are drawing a square
      predictions[75:80, 250:350, 300:400, 2] += create_circular_mask(100, 100,
      ↪center=None, radius=50).astype('int')
      predictions[75:80, 100:200, 200:300, 3] += create_circular_mask(100, 100,
      ↪center=None, radius=50).astype('int')
      predictions[75:80, 100:200, 200:300, 3] -= create_circular_mask(100, 100,
      ↪center=None, radius=33).astype('int')
      predictions[75:80, 100:200, 200:300, 3] += create_circular_mask(100, 100,
      ↪center=None, radius=15).astype('int')
```

```
[32]: display_slices(image, np.argmax(predictions, axis=-1), skip = 1) # visualize
      ↪our square on the image
```







Convert the NumPy arrays into RT-Structure

```
[33]: Dicom_reader.prediction_array_to_RT(prediction_array=predictions,
    ↪output_dir=RT_path,
    ROI_Names=['square', 'circle', 'target'])
```

```
Running off a template
Writing data for square
Writing data for circle
Writing data for target
Writing out data...Example_Data\RT_Structures
Finished!
```

2 Final notes

2.0.1 I hope that this code has been useful, if you have any suggestions or problems, please open an issue ticket or merge request on the Github: https://github.com/brianmanderson/Dicom_RT_and_Images_to_Mask

Thank you!

```
[ ]:
```